
Operating manual

Recirculating Coolers

FL1201 FL1203

FL1701 FLW1701 FL1703 FLW1703



Distributed by:



ADVANCED APPLIED TECHNOLOGIES Contact Us:

Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com

Changes without prior notification reserved Printed in Germany

1.951.4822 BE2 08/07 Druck: 15.08.2007

19514822_FL12xx_FL17xx.doc



JULABO Labortechnik GmbH 77960 Seelbach / Germany

+49 (0) 7823 / 51-0

+49 (0) 7823 / 24 91

info@JULABO.de

www.JULABO.de

Congratulations!

You have made an excellent choice.

JULABO thanks you for the trust you have placed in us.

This operating manual has been designed to help you gain an understanding of the principles of operating and possibilities of our circulators. For optimum utilization of all functions, we recommend that you thoroughly study this manual prior to beginning operation.

EC Conformity



The products described in the operating instructions conform to the requirements of the following European guidelines:

Low voltage regulations with respect to legal harmonization of the member countries concerning electric devices for use within certain voltage limits.

EMC guideline with respect to legal harmonization of the member countries concerning electromagnetic compatibility.



Quality Management System



The JULABO Quality Management System:

Development, production and distribution of temperature application instruments for research and industries conform to the requirements according to DIN EN ISO 9001:2000.

Certificate Registration No. 01 100044846

TABLE OF CONTENTS

| 1. Use | according to intended purpose | 4 |
|----------|------------------------------------------------------------|----|
| 2. Ope | rator responsibility – Safety recommendations | 4 |
| 2.1. | Handling | 4 |
| 2.2. | Use | 5 |
| 2.3. | Disposal | 5 |
| 3. Tech | hnical specifications | 6 |
| 3.1. | Warning functions and safety installations | 9 |
| 3.2. | Cooling water | 9 |
| 4. Safe | ety notes for the user | |
| 4.1. | Safety recommendations | 10 |
| • | acking and checking | |
| 6. Insta | allation | |
| 6.1. | • | |
| - | erating controls and functional elements | |
| • | rating procedures | |
| 8.1. | | |
| 8.2. | | |
| 8.3. | Filling | |
| 8.4. | Switching on / Start - Stop | |
| 8.5. | Setting the feed pressure | |
| 8.6. | Setting the temperatures | |
| 8.7. | AUTOSTART ON / OFF | |
| 8.8. | Remote control: activate – deactivate | |
| | ety installations | |
| 9.1. | | |
| 9.2. | • | |
| 10. | Troubleshooting guide / Error messages | |
| 11. | Electrical connections | |
| 12. | Remote control | |
| | . Setup for remote control | |
| | 2. Communication with a PC or a superordinated data system | |
| | 3. List of commands | |
| | Status messages | |
| | 5. Error messages | |
| 13. | Maintaining the cooling performance | |
| 14. | Adequate storing of operating manual | |
| 15. | Cleaning / repairing the unit | |
| | . Draining | |
| 16. | Warranty conditions | 27 |

1. Use according to intended purpose

JULABO recirculating coolers have been designed for temperature application to specific fluids. The pump connections can be used for cooling applications in an external circuit at a constant temperature.









- ☑ The recirculating coolers are operated via the splash-proof keypad. The implemented microprocessor technology allows to set and to store the setpoint that can be indicated on the LED temperature display.
- ☑ The PID temperature regulation is used to withdraw heat from the bath fluid by means of the cooling machine and to automatically regulate the required need.
- ☑ Electrical connections:
 - 1. The serial interface RS232 allows modern process technology without additional interface.
 - 2. Alarm output for external alarm message.
- Manually adjustable by-pass (handwheel) to reduce the pump capacity (e. g. for glass equipment).



JULABO recirculating coolers are not conceived for direct temperature application to food and luxury articles or pharmaceutical and medico-technical products. Direct temperature application means: Unprotected contact of the object with the bath medium (bath fluid).

2. Operator responsibility – Safety recommendations

The products of JULABO Labortechnik GmbH warrant a safe operation if installation, operation and maintenance is carried out according to common safety regulations. This section informs you about potential dangers that may arise from operating the recirculating cooler and also mentions the most important safety precautions.

The operator is responsible for the qualification of the personnel operating the units.

The operator should be constantly informed about the dangers involved with their job activities as well as preventive actions.

Make sure all persons expected to carry out operation, installation and maintenance of the unit read and understand the safety information and operating instructions.

When using hazardous materials, the circulator may only be operated by persons that are absolutely familiar with these materials and the circulator. These persons must be fully aware of possible risks. If you have any questions concerning the operation of your unit or the information in this manual, please contact us!

Contact JULABO

JULABO Labortechnik GmbH Eisenbahnstraße 45 77960 Seelbach / Germany

2 +49 7823 51-0

info@JULABO.de

2.1. Handling

You received a product conceived for industrial use. Nevertheless, avoid strikes to the housing, vibrations, damages to the keypad foil (keys, display) or contamination.

Make sure the product is regularly checked for proper condition. Regularly check (at least every 2 years) the proper condition of the mandatory, warning, prohibition and safety labels.

Take care that the mains supply features a low impedance to avoid any negative affects on the instrument being operated in the same mains.

This unit is designed for operation in a controlled electromagnetic environment. This means that transmitting devices (e.g. cellular phones) should not be used in the immediate vicinity. Magnetic radiation may influence other units with components susceptible to magnetic fields

(e.g. a monitor). We recommend to keep a minimum distance of 1 m.

Permissible ambient temperature: max. 40 °C, min. 5 °C.

Permissible relative air humidity: 50 % (40 °C).

Do not store in an aggressive atmosphere. Protect from contaminations. Do not expose to sunlight.

Appropriate Operation

Only qualified personnel is authorized to perform configuration, installation, maintenance and repairs of the recirculating cooler. Untrained personnel should be instructed by trained personnel.

2.2. Use

For the use according to the intended purpose, special material requirements have to be respected (bath fluids). Only use non-acid and non corroding materials. Observe all warnings for the used materials (bath fluids) and the respective instructions (safety data sheets).

Only use the unit in well ventilated areas. (see page 10).

The recirculating coolers are not for use in explosive atmosphere

When using hazardous materials, **the user must** attach the enclosed safety labels to the front of the unit so they are well visible.

Warning label W09:

Colors:

yellow, black

Mandatory label M018

Colors: blue, white

Semi S1-0701 Table A1-2 #9



Danger area.

Attention! Observe instructions.

(operating manual, safety data sheet)

Carefully read the user information prior to beginning operation

Scope: EU



Carefully read the user information prior to beginning

operation Scope: NAFTA

Observe the instructions in the manuals for instruments of a different make that you connect to the circulator, particularly the respective safety recommendations. Also observe the pin assignment of plugs and technical specifications of the products.

2.3. Disposal

This unit contains the refrigerants R134a or R404A – at this time considered not to have any negative effects on the ozone layer. However, during the long operating period of the unit, disposal prescriptions may change. So only qualified personnel should take care of disposal.



Directive 2002/96/EC of the European Parliament and of the Council of 27 January 2003 on waste electrical and electronic equipment (WEEE).

This directive requires electrical and electronic equipment marked with a crossedout trash can to be disposed of separately in an environmentally friendly manner.

Contact an authorized waste management company in your country.

Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!



3. Technical specifications

| Recirculating Cooler | | | FL | 1201 | FL1203 |
|----------------------------------|------------------|--------|-----|---------------------|------------------|
| Cooling | | | air | cooled | air cooled |
| Working temperature range | • | °C | -20 | +40 | -20 +4 0 |
| Temperature stability | | °C | | ±0.5 | |
| Temperature selection: | | | | digital | |
| via key pad | | | | indication on LEI | D-DISPLAY |
| remote control via pers | onal computer | | | indication on mo | nitor |
| Temperature indication: | | | | LED-DISPLAY | |
| Resolution | | °C | | 0.1 | |
| Temperature control | | | | PID 1 | |
| Temperature sensor | | | | Pt 100 | |
| Excess temperature protect | tion | | | 85 °C - fixed valu | ıe |
| Low liquid level protection | | | | float switch | |
| Cooling capacity | | °C | +20 | 0 -10 | <u>+20 0 -10</u> |
| Medium: Mixture water-glyd | ol | kW | 1.2 | 0.9 0.6 | 1.2 0.8 0.5 |
| Cooling compressor | | | | 1-stage | |
| Refrigerant | | | | R404A | |
| Electrical connections: | | | | | |
| Computer interface | | | | RS232 | |
| Alarm output | | | | for external alarn | n signal |
| · | | | | | |
| Circulating pump: | | | | | |
| discharge, max. | at 0 bar | I/min | 40 | | 60 |
| pressure, max. | at 0 liter | bar | 1.0 | | |
| pressure, adjustable | at 0 Liter | bar | | | 0.5 3.0 |
| Feed pressure indication | | bar | Ма | nometer | Manometer |
| Filling level indicator | | | | sight glass | |
| Filling volume from to | | liters | | 12 17 | |
| Dimensions (WxLxH) | | cm | | 50x76x64 | |
| Weight | | kg | 69 | | 73 |
| Ambient temperature range | • | °C | | 5 40 | |
| Return flow temperature | | °C | | 80 max. | |
| IP class according to IEC 60 529 | | | | IP 21 | |
| Mains power connection 23 | 0 V/50 Hz | V/ Hz | 207 | '-253 / 50 | 207-253 / 50 |
| Current input at 230 V | | Α | 6 | | 9 |
| Mains power connection 2 | 08 - 220 V/60 Hz | V/ Hz | | ' - 242 / 60 | 197 - 242 / 60 |
| Current input at 208V / 220 | | A | | | |
| Mains power connection 1 | 15 V/60 Hz | V/ Hz | 103 | 3 -127 / 60 | 103 -127 / 60 |
| Current input at 115 V | | Α | | | |
| · | | | | | |

All measurements have been carried out at: rated voltage and frequency, ambient temperature: 20 $^{\circ}\text{C}$

| Recirculating Cooler | | | FL | 1701 | FLW1701 |
|----------------------------------|---------------------------|--------|-----|--------------------|-----------------|
| Cooling | | | air | cooled | water cooled |
| Working temperature range | | °C | -20 | +40 | -20 +4 0 |
| Temperature stability | | °C | | ±0.5 | |
| Temperature selection: | | | | digital | |
| via key pad | | | | indication on LE | ED-DISPLAY |
| remote control via pers | onal computer | | | indication on mo | onitor |
| Temperature indication: | | | | LED-DISPLAY | |
| Resolution | | °C | | 0.1 | |
| Temperature control | | | | PID 1 | |
| Temperature sensor | | | | Pt 100 | |
| Excess temperature protect | tion | | | 85 °C - fixed va | lue |
| Low liquid level protection | | | | float switch | |
| Cooling capacity | | °C | +20 | 0 0 -10 | +20 0 -10 |
| Medium: Mixture water-glyd | ol | kW | 1.7 | 1.1 0.85 | 1.7 1.1 0.85 |
| Cooling compressor | | | | 1- stage | |
| Refrigerant | | | | R404A | |
| | | | | | |
| Electrical connections: | | | | | |
| Computer interface | | | | RS232 | |
| Alarm output | | | | for external alarm | signal |
| Circulating pump: | | | | | |
| discharge, max. | at 0 bar | I/min | 40 | | 40 |
| pressure, max. | at 0 liter | bar | 1.0 | | 1.0 |
| pressure, adjustable | at 0 Liter | | | | |
| Feed pressure indication | | bar | Ма | nometer | Manometer |
| Filling level indicator | | | | sight glass | |
| Filling volume from to | | liters | | 12 17 | |
| Dimensions (WxLxH) | | cm | | 50x76x64 | |
| Weight | | kg | 71 | | 72 |
| Ambient temperature range | ı | °C | | 5 40 | |
| Return flow temperature | | °C | | 80 max. | |
| Cooling water | | | | | |
| Flow rate at 20 °C inlet | temperature | I/min | | - | 2.8 |
| IP class according to IEC 60 529 | | | | IP 21 | |
| | 0.1450 | | | 7 050 / 50 | 007.050.450 |
| Mains power connection 23 | U V/50 Hz | V/ Hz | _ | 7-253 / 50 | 207-253 / 50 |
| Current input at 230 V | 10 220 1 <i>11</i> 60 LI- | Α | 6 | - 040400 | 6 |
| Mains power connection 20 | | V/ Hz | 197 | 7 - 242 / 60 | 197 - 242 / 60 |
| Current input at 208V / 220 | | A | 400 | . 407 / 00 | 400 407 / 00 |
| Mains power connection 1 | IO V/OU FIZ | V/ Hz | 103 | 3 -127 / 60 | 103 -127 / 60 |
| Current input at 115 V | | Α | | | |

All measurements have been carried out at: rated voltage and frequency, ambient temperature: 20 $^{\circ}\text{C}$

| Recirculating Cooler | | | FL1703 | FLW1703 |
|----------------------------------------------------|-----------------|--------|------------------|------------------|
| Cooling | | ••• | air cooled | water cooled |
| Working temperature range Temperature stability | | °C | -20 +40 | -20 + 40 |
| • | | °C | ±0.5 | |
| Temperature selection: | | | digital | |
| via key pad | | | | LED-DISPLAY |
| remote control via perso | onal computer | | indication or | n monitor |
| Temperature indication: | | | LED-DISPLA | AY |
| Resolution | | °C | 0.1 | |
| Temperature control | | | PID 1 | |
| Temperature sensor | | | Pt 100 | |
| Excess temperature protect | ion | | 85 °C - fixed | value |
| Low liquid level protection | | | float switch | |
| Cooling capacity | | °C | +20 0 -10 | <u>+20 0 -10</u> |
| Medium: Mixture water-glyc | ol | kW | 1.7 1.0 0.75 | 1.7 1.0 0.75 |
| Cooling compressor | | | 1- stage | |
| Refrigerant | | | R404A | |
| Electrical connections: | | | | |
| Computer interface | | | RS232 | |
| Alarm output | | | for external ala | rm signal |
| Circulating pump: | | | | |
| discharge, max. | at 0 bar | I/min | 60 | 60 |
| pressure, max. | at 0 liter | bar | 0.5 3.0 | 0.5 3.0 |
| pressure, adjustable | at 0 Liter | | | |
| Feed pressure indication | | bar | Manometer | Manometer |
| Filling level indicator | | | sight glass | |
| Filling volume from to | | liters | 12 77 | |
| Dimensions (WxLxH) | | cm | 50x76x64 | |
| Weight | | kg | 76 | 77 |
| Ambient temperature range | | °Č | 5 40 | |
| Return flow temperature | | °C | 80 max. | |
| Cooling water | | | | |
| Flow rate at 20 °C inlet | temperature | I/min | | 2.8 |
| IP class according to IEC 60 | | | IP 21 | |
| Mains power connection 236 | 0 V/50 Hz | V/ Hz | 207-253 / 50 | 207-253 / 50 |
| Current input at 230 V | | Α | 9.0 | 9.0 |
| Mains power connection 20 | 8 - 220 V/60 Hz | V/ Hz | 197 - 242 / 60 | 197 - 242 / 60 |
| Current input at 208V / 220 | | Α | | |
| Mains power connection 11 | | V/ Hz | 103 -127 / 60 | 103 -127 / 60 |
| Current input at 115 V | | Α | | |

All measurements have been carried out at: rated voltage and frequency, ambient temperature: 20 $^{\circ}\text{C}$

3.1. Warning functions and safety installations

85 °C - fixed value Excess temperature protection

Low liquid level protection float switch

Alarm message optical + audible (permanent)

Excess temperature - Warning function 75 °C

Overload protection for compressor and pump motor

Classification according to DIN 12876-1 class I

Environmental conditions according to IEC 61 010-1:

Use only indoor.

Altitude up to 2000 m - normal zero.

Ambient temperature: +5 ... +40 °C (for storage and transportation)

Air humidity:

Max. rel. humidity 80 % for temperatures up to +31 °C,

linear decrease down to 50 % relative humidity at a temperature of +40 °C

Power supply: corresponds to Class I; according to VDE 0106 T1

2



Not for use in explosive atmosphere

Overvoltage category Ш Pollution degree

3.2. **Cooling water**

Only for water cooled models - FLW:

Cooling water pressure (IN / OUT) 4,5 bar max.

Difference pressure (IN - OUT) 2 bar bis 4,5 bar

Flow rate 2,8 I/min typical Cooling water temperature 20 °C

Quality of cooling water:

pH at 25 °C 7 to 8.5 <30 mg/l Suspended matter Size of suspended matter 0.1 mm max.

Growth of algae not permissible

4. Safety notes for the user



In addition to the safety warnings listed above, warnings are posted throughout the manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)."

The danger is described according to an alarm keyword.

Read and follow these important instructions.



Warning:

Describes a possibly highly dangerous situation. If this is not avoided, serious injury and danger to life could result.



Caution:

Describes a possibly dangerous situation. If this is not avoided, slight or minor injuries could result.

A warning of possible damage can also be contained in the text.



Notice:

Describes a possibly harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.

4.1. Safety recommendations

Follow the safety recommendations to prevent damage to persons or property. Further, the valid safety instructions for working places must be followed.



- Only connect the unit to a power socket with earthing contact (PE protective earth)!
- Place the instrument on an even surface on a pad made of non-inflammable material.
- Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your unit.
- Never operate the unit without bath fluid in the bath.
- Exercise caution when emptying hot bath fluids!
 Check the temperature of the bath fluid prior to draining (by switching the unit on for a short moment for example).
- Employ suitable connecting tubing.
 - Make sure that the tubes are securely attached.
- Never operate damaged or leaking equipment.
- Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures, or before moving the unit.



- Always empty the bath before moving the unit.
- Never operate equipment with damaged mains power cables.
- Risk of injury for hands. Close cover carefully.

5. Unpacking and checking

Unpack the recirculating cooler and accessories and check for damages incurred during transit. These should be reported to the responsible carrier, railway, or postal authority, and a request for a damage report should be made. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.




← Lifting device for transportation by crane see page 13


6. Installation

Rear view:

- Place the unit on an even surface on a pad made of non-flammable material.
- Cooling machine, pump motor and electronics produce intrinsic heat that is dissipated via the venting openings.! Never cover these openings!
- Keep at least 20 cm of open space on the front and rear venting grids.
 - Do not set up the unit in the immediate vicinity of heat sources and do not expose to sun light.



Example: FL1201 (air cooled)




The place of installation should be large enough and provide sufficient air ventilation to ensure the room does not warm up excessively because of the heat the instrument radiates to the environment. (Max. permissible ambient temperature: 40 °C). With regard to a disturbance in the cooling loop (leakage), the guideline EN 378 prescribes a certain room space to be available for each kg of refrigerant. > For 0.48 kg of refrigerant R404A, a room space of 1 m³ is required.

Example: model FL1703with 0.7 kg filling quantity of refrigerant R404A = 1.46 m³ volume


Connect the tubings for cooling the external system to the pump connectors M16x1 for feed and return (14) on the rear of the recirculating cooler.



- Connect a piece of tubing to the overflow connector (15) and drain into a suitable vessel, which always has to be placed lower thant the exit "Overflow".
- Turn the adjusting valve (14) counter-clockwise to set the lowest manometric pressure.
- Before operating the unit after transport, wait about one hour after setting it up. This will allow any oil that has accumulated laterally during transport to flow back down thus ensuring maximum cooling performance of the compressor.



Example: FL1703 (watter cooled, pressure, adjustable 0,5 ... 3, 0 bar)



Only water cooled models - FLW: Ensure circulation of cooling water by connecting the tubing to cooling water inlet (IN)and outlet (OUT) on the rear (16) of the recirculating cooler.

Cooling water see page 9.

G3/4" external thread Cooling water connectors 12 mm inner dia. tubing **Tubing**

Cooling water inlet OUT Cooling water outlet



Notice: Cooling water circuit

Risk of oil leaking from the cooling circuit (compressor) of the recirculating cooler into the cooling water in case of a fault in the circuit!

Observe the laws and regulations of the water distribution company valid in the location where the unit is operated.



Caution: Pump pressure

- Determine and check the max. admissible pressure for the external circuit before
 putting the unit into operation. The max. pressure is determined by the weakest
 element in the circuit (e. g. glass equipment).
- · Securely attach all tubing to prevent slipping.



Notice: Flood hazard!.

In case the system to be cooled is located at a higher level than the recirculating cooler, take note of bath liquid flowing back when the unit is switched off.

Return flow safety device

Should the filling volume of the bath tank not be sufficient, prevent the liquid from flowing back by using shut-off valves..

| Order No. | Description | Suitable for |
|-----------|----------------------------------------|--------------------|
| 8 970 456 | Shut-off valve for loop circuit, M16x1 | FL1201 / FL(W)1701 |
| 8 970 454 | Shut-off valve G ¾" | FL1203 / FL(W)1703 |

The following questions shall help to recognize possible dangers and to reduce the risks to a minimum.

- Are all tubes and electrical cables connected and installed?
 Note:
 - sharp edges, hot surfaces in operation, moving machine parts, etc.
- What to do when a dangerous substance was spilled on or in the unit?
 Before starting to work, obtain information concerning the substance and determine the method of decontamination.

6.1. Tubing




Caution:

- Employ suitable connecting tubing.
- Make sure that the tubing is securely attached.
- Avoid sharp bends in the tubing, and maintain a sufficient distance from surrounding walls.
- Regularly check the tubing for material defects (e.g. for cracks).
- Preventive maintenance: Replace the tubing from time to time.


Recommended tubing:

| Order No. | Description | Suitable for |
|--------------|----------------------------------------------------------|------------------------------------|
| 8930308 | 1 m CR [®] -tubing 8 mm inner dia. (-20 +120°C) | FL1201 |
| 8930312 | 1 m Reinforced tubing 8 mm inner dia. (-40 +12 | 20°C) FL1201 |
| 8930319 | 1 m Reinforced tubing 12 mm inner dia. (-40 +12 | 20°C) FL(W)1203/1703 |
| Tubing insul | ation | |
| 8930412 | 1 m Insulation, 18 mm inner dia. | Reinforced tubing 8 mm inner dia. |
| 8930413 | 1 m Insulation, 23 mm inner dia. | Reinforced tubing 12 mm inner dia. |
| 8930419 | 1 m Insulation, 29 mm inner dia. | Reinforced tubing ¾" inner dia. |
| Tube clamps | • | |
| 8970481 | 2 Tube clamps, size 2 | Reinforced tubing 8 mm inner dia. |
| 8970482 | 2 Tube clamps, size 3 | Reinforced tubing 12 mm inner dia. |
| 8970483 | 2 Tube clamps, size 4 | Reinforced tubing ¾" inner dia. |
| | | |

7. Operating controls and functional elements



| 1 | | Mains power switch, spash-water protected on |
|-----|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|
| | 0 | O off |
| 2.0 | | Keypad, spash-water protected |
| 2.1 | | Edit keys (set point increase or decrease) |
| 2.2 | ↓ | Enter key Store set point value / parameter |
| 3.0 | | Indication |
| 3.1 | - 15.0 | LED temperature display |
| 3.2 | ₩. | Control indicator – Cooling |
| 3.3 | lack | Control indicator – Alarm |
| 4 | | Protection lid for fill in opening |
| 5 | ∘ (*****) ∘ RS232 | Interface RS232: remote control via personal computer |
| 6 | | Alarm output (for external alarm signal) |
| 7 | | Feed pressure indication: Manometer |
| 8 | | Filling level indication |
| 9 | | Drain tap with drain port |
| 10 | | Venting grid, removable |
| 11 | 16 M M | Mains circuit breakers (resettable) |
| 12 | | Mains power cable with plug |
| 13 | IN OUT | Pump connector OUT – pressure pump |
| | 0 0 | Pump connector IN– Return FL/FLW1201 M16x1 FL/FLW1203 G ¾ external thread FL/FLW1703 G ¾ external thread |
| 14 | | Feed pressure adjusting valve |
| 15 | | Overflow connector |
| 16 | ST CONTRACTOR OF THE PROPERTY | Only water cooled models IN Cooling water inlet OUT Cooling water outlet G3/4" external thread 12 mm inner dia. tubing |
| 17 | ~ | Castor (at the back) |
| 18 | | Protection lid for storing place of operating manual |



- 4 Protection lid for fill in opening
- 18 Protection lid for storing place of operating manual

8. Operating procedures

8.1. Bath fluids



Caution:

No liability for use of other bath liquids!

Do not use alcohols.

Water:

The quality of water depends on local conditions. Ferrous water can cause corrosion - even on stainless steel. Chloric water can cause pitting corrosion.

Water: - No liablity for use with water.

Danger of freezing at working temperatures <5 °C.

Mixture water -glycol:

Strictly observe the safety data and handling instructions from the manufacturer.

The proportion of water might evaporate by and by. Check the mixing ratio regularly and refill water if necessary.

Recommended bath fluids:

| Bath fluids | Temperature range |
|------------------------------|-------------------|
| Julabo Thermal G | -30 °C 80 °C |
| mixture water/glycol (50:50) | -30 °C 50 °C |
| water | +5 °C 80 °C |

| Order No. | Desription | Quantity |
|-----------|------------------|-----------|
| 8 940 124 | JULABO Thermal G | 10 liters |
| 8 940 125 | JULABO Thermal G | 5 liters |



Notice:

Please contact JULABO before using other than recommended bath fluids. JULABO takes no responsibility for damages caused by the selection of an unsuitable bath fluid

8.2. Power connection



Caution:

Only connect the unit to a power socket with earthing contact (PE – protective earth)! We disclaim all liability for damage caused by incorrect line voltages!

Check to make sure that the line voltage matches the supply voltage specified on the identification plate.


8.3. Filling



Notice:

Risk of injury for hands. Close cover carefully.

Top view



Take care that no liquid enters the interior of the circulating cooler.


(i) Connect the tubing from the external system to the pump connectors and check for leaks



Respect instructions from page 10 to page 12!

- (i) Check to make sure that the drain tap (9) is closed.
- Unlock and open lid of fill in opening (4) by slightly pushing.
- Fill in tempering fluid up to marking "H" of the filling level indicator.
- Turn the mains switch (1) on (Switching on see page 17)
- Switch on unit. To do so press button for approx. 4 seconds.
- Tempering fluid is pumped into the externally connected system. Refill fluid up to marking "H".
- The recirculating cooler is ready for operation.

8.4. Switching on / Start - Stop



Switching on:

The recirculating cooler is turned on and off with the mains switch. The unit performs a self-test. All segments of the 4-digit LED temperature DISPLAY and all indicator lights will illuminate (as illustrated on the left).

Then the software version and the type of unit is indicated. Examples: (v 1.02) (FL1703)

The display "**OFF**" indicates the unit is ready to operate (standby mode).



Start: Press enter for about 4 seconds.

The LED temperature DISPLAY indicates the actual bath

temperature.

Stop: Press enter for about 4 seconds.

Turn the unit off with the mains power switch.

8.5. Setting the feed pressure





Set the max. permissible feed pressure (example: 2 bar) by slowly closing the adjusting valve (14) on the rear and looking at the manometer (7).

The max. pressure is determined by the weakest element in the circuit (e. g. glass equipment).

8.6. Setting the temperatures

- Press one of the keys for a short moment.
 The setpoint value instead of the actual value is indicated on the display for about 8 seconds.
 The value can now be changed.
- 2. Change value:

Press <u>to set a higher value.</u>

Press to set a lower value.

Keep the keys depressed for the value to change fast.

3. Press enter to store the value.

8.7. AUTOSTART ON / OFF

The recirculating cooler has been configured and supplied by JULABO according to N.A.M.U.R. recommendations. This means for the start mode, that the unit must enter a safe operating state after a power failure (non-automatic start mode). This safe operating state is indicated by "OFF" on the LED temperature display. A complete shutdown of the main functional elements such as compressor and circulating pump is effected simultaneously.

Should such a safety standard not be required, the AUTOSTART function (automatic start mode) may be activated, thus allowing the start of the circulator directly by pressing the mains power switch or using a timer.



Keep depressed enter and turn on the unit with the mains power switch.

For a short while the LED DISPLAY indicates the effective start mode:








Warning:

For supervised or unsupervised operation with the AUTOSTART function, avoid any hazardous situation to persons or property.

The circulator does no longer conform to N.A.M.U.R. recommendations.

8.8. Remote control: activate – deactivate



(Interface OFF)





The recirculating cooler is to be prepared for remote control by a personal computer via the serial interface RS232. Set the interface item from >IOFF< to >ION<.

Remote control: activate - deactivate:

- Switch off recirculating cooler by pressing the mains switch and wait approx. 5 seconds.
- Keep depressed the keys and enter simultaneously and turn on the unit with the mains power switch.
- >I OFF< No remote control via RS232 (Factory setting)
- >I On< Remote control via RS232
- (i) The software version and the type of unit is indicated (see example on the left).

The display "rOFF" indicates the unit is ready to be operated via remote control.

9. Safety installations

9.1. Excess temperature protection



This safety installation is independent of the control circuit.

When the temperature of the bath fluid has reached the safety temperature (85 °C), a complete shutdown of the compressor and pump is effected.

The alarm is indicated by optical and audible signals (continuous tone) and on the LED-DISPLAY appears the error message "Error 14".

9.2. Low level protection



This safety installation is independent of the control circuit.

If the low liquid level protection device is triggered, a complete shutdown of the compressor and circulating pump is effected.

The alarm is indicated by optical and audible signals (continuous tone) and on the LED-DISPLAY appears the error message "Error 01".

Turn off the unit with the mains switch, refill bath fluid and turn the unit on again!



Caution:

For refill always use the same bath fluid type that is already in the bath.



Notice:

Check the low liquid level protection device at least twice a year!

To execute a functional test, drain the liquid until the alarm for low liquid level is triggered. Refill liquid afterwards.

10. Troubleshooting guide / Error messages



Whenever the microprocessor electronics registers a failure, a complete shutdown of the compressor and circulating pump is performed. The alarm light

"A" illuminates and a continuous signal tone sounds.

The LED temperature display indicates the cause for the alarm in form of a code.



Press enter to quit the audible signal.

- The recirculating cooler is operated without bath fluid, or the liquid level is insufficient.
 - Replenish the bath tank with the bath fluid.
- Tube breakage has occured (insufficient filling level due to excessive bath fluid pumped out). Replace the tubing and replenish the bath tank with the bath fluid.



Cable of the working temperature sensor interrupted or short-circuited.

| E 12 | Error in A/D converter |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| E 14 | The return temperature is above the switch-off value of the high temperature protection (85°C). Check dimensioning of application. Use a stronger recirculating cooler if necessary. |
| | After eliminating the malfunction, press the mains power switch off and on again to cancel the alarm state. If the unit cannot be returned to operation, contact an authorized service station. |
| E 03 | Warning without a complete shutdown of the unit Excess temperature warning starting at 75 °C The return temperature soon reaches the swith-off value of the high temperature protection (85 °C). |
| E 20 | Cooling of the condenser is affected. (see page 24) Clean air-cooled condenser. Check the flow rate and cooling water temperature on water-cooled condenser. |

If the unit cannot be returned to operation, contact an authorized JULABO service station.

Disturbances that are not indicated.

Overload protection:: a) for cooling machine

b) for pump motor

After a short cooling interval, the unit will automatically start running.

16 A

Mains circuit breakers (resettable) –10A.

11. Electrical connections



Notice: Use shielded cables only.

The shield of the connecting cable is electrically connected to the plug housing.

The unit ensures safe operation if connecting cables with a maximum length of 3 m are used. The use of longer cables does not affect proper performance of the unit, however external interferences may have a negative impact on safe operation.

RS232 serial interface

This port can be used to connect a computer with an RS232 cable for remote control of the recirculating cooler.


RS232C

Pin assignments:

| in 2 | RxD | Receive Data |
|-------|------|---------------------|
| Pin 3 | TxD | Transmit Data |
| Pin 5 | 0 VD | Signal GND |
| Pin 6 | DTR | Data terminal ready |
| Pin 7 | RTS | Request to send |
| Pin 8 | CTS | Clear to send |

Accessories:

| Order No. | Description |
|-----------|---------------------------------------------|
| 8 980 073 | RS232 interface cable 9-pol./9-pol. , 2,5 m |
| 8 900 110 | USB interface adapter cable |



Alarm output


Potential-free change-over contact for external alarm signal.

Pin 2 and 3 are connected in case of an alarm. Pin 2 and 1 are connected in "OFF" or "rOFF" condition or mains switch "Off".

Switching capacity max. 30 W / 40 VA max. 125 V~/-Switching voltage Switching current max. 1 A

12. Remote control

12.1. Setup for remote control



Check the interface parameters for both interfaces (on recirculating cooler and PC) and make sure they match.

Interface parameters are pre-determined.

Type RS232 Baudrate 4800 bauds **Parity** even

Handshake hardware handshake

12.2. Communication with a PC or a superordinated data system

If the recirculating cooler is put into remote control mode the MULTI-DISPLAY (LED) will read "R -OFF-" = REMOTE STOP. The recirculating cooler is now operated via the computer.

In general, the computer (master) sends commands to the recirculating cooler (slave). The recirculating cooler sends data (including error messages) only when the computer sends a query.

In remote control mode:

After a power interruption the order to start and all values which have to be adjusted must be resent from the personal computer via the interface. AUTOSTART is not possible.

A transfer sequence consists of:

command out/in command

space (⇔; Hex: 20) out/in command

• parameter (the character separating decimals in a group is the

period) out command

end of file (∠; Hex: 0D) out/in command

 The response (data string) after an in command is always followed by a line feed (LF, Hex: 0A).

Important times for a command transmission:

To ensure a safe data transfer, the time gap between two commands should be at least 250 ms.

The recirculating cooler automatically responds to an **in** command with a data string followed by a LF (Line Feed). The next command should only be sent after 10 ms.

The commands are divided into in or out commands.

in commands: asking for parameters to be displayed

out commands: setting parameters

The **out** commands are valid only in remote control mode.

Examples:

Command to set the working temperature to 15,5 °C:

out_sp_00 ⇔ 15.5↓

Command to ask for the working temperature

in sp 00↓

Response from the recirculating cooler:

15.5↓ LF

(S)

12.3. List of commands

out commands: Setting parameters or temperature values.

| Command | Parameter | Response of recirculating cooler |
|-------------|-----------|----------------------------------|
| out_mode_05 | 0 | Stop the unit = R –OFF |
| out_mode_05 | 1 | Start the unit. |
| out_sp_00 | xxx.xx | Set working temperature |

in commands: Asking for parameters or temperature values to be displayed.

| Command | Parameter | Response of recirculating cooler |
|------------|-----------|------------------------------------------------------------------|
| version | none | Number of software version (V X.xx) |
| status | none | Status message, error message (see page 23) |
| in_pv_00 | none | Actual bath temperature. |
| in_sp_00 | none | Working temperature |
| in_mode_05 | none | Recirculating cooler in Stop/Start condition: 0 = Stop 1 = Start |

12.4. Status messages

| Status messages | Description |
|-----------------|----------------------------------------------|
| 00 MANUAL STOP | Recirculating cooler in "OFF" state. |
| 01 MANUAL START | Recirculating cooler in keypad control mode. |
| 02 REMOTE STOP | Recirculating cooler in "r OFF" state. |
| 03 REMOTE START | Recirculating cooler in remote control mode. |

12.5. Error messages

| Error messages | Description |
|---------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| -01 LOW LEVEL ALARM | Low liquid level alarm. |
| -05 WORKING SENSOR ALARM | Working temperature sensor short-circuited or interrupted. |
| -03 EXCESS TEMPERATURE WARNING | High temperature warning. Starting at 75 °C (no deactivation) The return temperature soon reaches the switch-off value of the high temperature warning function (85 °C) |
| -07 I ² C-BUS ERROR | Internal error when reading or writing the I ² C bus. |
| -08 INVALID COMMAND | Invalid command. |
| -09 COMMAND NOT ALLOWED IN CURRENT OPERATING MODE | Invalid command in current operating mode. |
| -10 VALUE TOO SMALL | Entered value too small. |
| -11 VALUE TOO LARGE | Entered value too large. |


| Error messages | Description |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| -12 TEMPERATURE MEASUREMENT ALARM | Error in A/D converter. |
| -14 EXCESS TEMPERATURE PROTECTOR ALARM | The return temperature is above the switch-off value of the high temperature warning function of 85 °C. Check dimensioning of application. Use a stronger recirculating cooler if necessary. |
| -20 WARNING: CLEAN CONDENSOR OR CHECK COOLING WATER CIRCUIT OF REFRIGERATOR | Cooling of the condenser is affected. Clean air-cooled condenser. Check the flow rate and cooling water temperature on water-cooled condenser. |

13. Maintaining the cooling performance



Notice:


Risk of injury for hands when mounting the venting grid.



Air cooled models = FL

To maintain the full cooling performance, clean the condenser from time to time.

- Switch off the unit, disconnect mains power cable.
- Hold the venting grid, pull out and remove.
- Clean the ribbed condenser with a vacuum cleaner.
- Replace the venting grid.
- Switch on the unit.



Water cooled models = FLW

In order to maintain a good condition of the cooling compressor, the sieve in the cooling water input should be cleaned in regular intervals.

- Switch the unit off, disconnect the power plug.
- Interrupt the cooling water input.
- Disconnect the tubing from the nozzle "IN" and take out the dirty sieve.
- Clean the sieve.
- Put in the sieve and reconnect the tubing.
- · Open the cooling water input.
- Take care the tubing connection is not leaking.
- The unit is ready to operate again.

14. Adequate storing of operating manual

Store the operating manual at the foreseen place at the unit and lock it by means of the protection lid (18).

15. Cleaning / repairing the unit



Caution:

Always turn off the unit and disconnect the mains cable from the power source before cleaning the unit.

Prevent humidity from entering into the circulator.

Electrical connections and any other work must be performed by qualified personnel only.

Cleaning:

Clean the outside of the unit using a wet cloth and low surface tension water.

The recirculating cooler is designed for continuous operation under normal conditions. Periodic maintenance is not required.

The tank should be filled only with a bath fluid recommended by JULABO. To avoid contamination, it is essential to change the bath fluid from time to time.

Repairs:

Before asking for a service technician or returning a JULABO instrument for repair, please contact an authorized JULABO service station.

When returning the unit:

- Clean the unit in order to avoid any harm to the service personnel
- Attach a short fault description.
 If you intend to return your JULABO unit to us, you will find a Service Return Form on our website www.julabo.de. Please use this as a delivery note and include it to the unit or send it in advance either by Fax or E-Mail.
- When returning a unit, take care of careful and adequate packing.
- JULABO is not responsible for damages that might occur from insufficient packing.



JULABO reserves the right to carry out technical modifications with repairs for providing improved performance of a unit.

15.1. Draining



Notice:

Store and dispose the used bath fluid according to the laws for environmental protection.



Risk of injury for hands when mounting the venting grid.



- Turn off the unit and disconnect the mains cable from the power source.
- Hold the venting grid, pull out and remove.
- Slide a short piece of tube onto the drain port and hold it into a pail.
- Open the drain tap and empty the unit completely.
- Close the drain tap and replace the venting grid.

16. Warranty conditions

JULABO Labortechnik GmbH warrants its products against defects in material or in workmanship, when used under appropriate conditions and in accordance with appropriate operating instructions

for a period of ONE YEAR.

Extension of the warranty period – free of charge



With the '1PLUS warranty' the user receives a free of charge extension to the warranty of up to 24 months, limited to a maximum of 10 000 working hours.

To apply for this extended warranty the user must register the unit on the JULABO web site www.julabo.de, indicating the serial no. The extended warranty will apply from the date of JULABO Labortechnik GmbH's original invoice.

JULABO Labortechnik GmbH reserves the right to decide the validity of any warranty claim. In case of faults arising either due to faulty materials or workmanship, parts will be repaired or replaced free of charge, or a new replacement unit will be supplied.

Any other compensation claims are excluded from this guarantee.