

Data sheet

Bst DNA Polymerase (Exonuclease Minus)

Cat. No: P0045 (2000 U)
Cat. No: P0046 (10,000 U)

Introduction:

Bst DNA Polymerase, Exonuclease Minus, is a 67 kDa *Bacillus stearothermophilus* DNA Polymerase protein (large fragment) which has a 5'-3' polymerase activity and strand displacement activity but lacks 3' – 5' exonuclease activity. Also has reverse transcription activity.

Source: A recombinant *E. coli* strain carrying the Bst DNA Polymerase gene (large fragment).

Application:

- ✓ nucleic acid amplification methods, including isothermal amplification
- ✓ whole genome amplification
- ✓ multiple displacement amplification
- ✓ sequencing DNA with high GC content and secondary structures
- ✓ rapid sequencing from nanogram amounts of DNA Template

Kit Contents

	P0045	P0046
Bst DNA Polymerase (8 U/μL)	250μl	5X250μl
Reaction Buffer (10x)	1.2 mL	5X1.2 mL
100 mM MgSO4	500 μl	5X500 μl

Storage:

Store at -20°C.

Unit definition:

1 unit is defined as the amount of polymerase required to convert 10 nmol of dNTPs into acid insoluble material in 30 minutes at 65°C.

Quality control:

- ✓ **Bst DNA Polymerase** is free of detectable RNase, and DNase (exo- and endonuclease) activities.
- ✓ Purity: >99% as judged by SDS-polyacrylamide gels with blue staining.
- ✓ No detectable DNA contamination

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk

 **canvax**

Typical LAMP Protocol

1.- Prepare the reaction mix as shown in the Table in the order listed. During this step the reaction mix tube should always be held on the ice to prevent the background activity of the enzyme.

Component	Volume (μ L)	Final Conc
10X ThermoPol Buffer	2.5 μ l	1X (contains 2 mM MgSO ₄)
MgSO ₄ (100 mM)	1.5 μ l	6 mM (8 mM total)
dNTP Mix (10 mM)	3.5 μ l	1.4 mM each
FIP/BIP Primers	-	1.6 μ M each
F3/B3 Primers	-	0.2 μ M each
Loop F/B Primers	-	0.4-0.8 μ M each
Bst DNA Polymerase 8 U/ μ l)	1 μ l	0.32 U/ μ l
DNA Sample	variable	> 10 copies or more
Nuclease-free Water	to 25 μ l	

2.-Incubate reaction at 65°C for 30–60 minutes. Running a temperature gradient from (55–65 °C) is strongly recommended to determine optimum temperature.

- Running a no-template control is strongly recommended to ensure amplification specificity.
- If optimization is desired, try titrating Mg²⁺ (4–10 mM final) or Bst DNA Polymerase (0.04–0.32 U/ μ l), or changing reaction temperature (50–68°C).

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for Material Safety Data Sheet of the product.

Some uses for this product may require licenses. Canvax Biotech does not encourage or support the unauthorized or unlicensed use of patented nucleic acid amplification processes for isothermal amplification, whole genome amplification (WGA), multiple displacement amplification (MDA), and Next Generation sequencing. It is the sole responsibility of the buyer to ensure that use of the product does not infringe the patent rights of third parties.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk