

Data sheet

Clean-Easy PCR Purification Kit

Cat. No: AN0063-S (20 reactions)

Cat. No: AN0063 (50 reactions)

Cat. No: AN0064 (100 reactions)

Description

Clean-Easy PCR Purification kit provides a rapid and efficient method to purify DNA and remove contaminants from reaction mixtures (e.g. PCR, digestion or labeling reactions). Clean-Easy minispin columns contain an exclusive membrane that allows DNA adsorption in presence of chaotropic salts and the removal of contaminants.

Features

- **Simple** and **Just a few minutes** procedure.
- **70-90% DNA recovery.**
- Suitable for DNA fragments as short as 75 bp.
- DNA purified **Ready to use** in all molecular biology procedures.

Applications

- Removal of proteins and salts from PCR, restriction digestion, dephosphorylation, ligation or labelling reactions.
- Changing of a restriction enzyme buffer.
- Re-purification of genomic DNA

Kit Components

ITEM	AN0063	AN0064
Clean-Easy minispin columns	50	100
Collection tubes (2 mL)	50	100
PB Buffer	25 ml	50 ml
PE Buffer *	11.25 ml	22.5 ml
EB Buffer	9 ml	9 ml

* Add the volume ethanol (96%-100%) specified [Not included] prior to initial use (see bottle label for volume). After ethanol has been added, mark the bottle to indicate that this step has been completed.



Store at RT.

Quality Certifications

Clean-Easy PCR Purification Kit is tested in the purification of a 0.6 kb DNA fragment from PCR mixture. The purified band is analysed in agarose gel electrophoresis.

(Continued on reverse side)

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432

Fax: (01) 4523967

E-mail: info@labunlimited.com

Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30

Fax: 08452 30 50 30


E-mail: info@labunlimited.co.uk

Web: www.labunlimited.co.uk

 **canvax**

Assay procedure

1. Add 5 volumes of **Buffer PB** to one volume of PCR solution and mix thoroughly by pipette.
2. Label the lid of a new spin column placed in a 2 ml collection tube. Carefully apply the mix from step 1 to the spin column and Centrifuge at 13000 rpm for 1 minute.
3. Place the spin column in a new 2 ml collection tube, and discard the collection tube containing the filtrate.
4. Add 700 µl of **buffer PE** for Wash to the minispin column and centrifuge at 13000 rpm for 1 minute.
Remember! Before using it for the first time, add ethanol (96–100%) to the PE Buffer as indicated on the bottle.
5. Discard the flow-through and centrifuge at 13000 rpm for 1 minute. *This step is essential for removing traces of PE buffer.*
6. Transfer the minispin column into a new, labeled 1.5 ml microcentrifuge tube.
7. Carefully open the minispin column and pipet 30 µl **Buffer EB** or H₂O (pH=7.0-8.5) directly onto the membrane. Close the cap and incubate for 1 min at room temperature, then centrifuge at 13000 rpm for 1 min to elute DNA. *To increase the DNA yield you can warm the buffer EB/H₂O to 65 °C before adding to the column.*



PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, and is not suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for the Material Safety Data Sheet of the product.

Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk