Nitrogen Generators

for Gas Chromatography

Distributed by:



Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com



Nitrogen on Demand, up to 2,000 ml/min

The Parker Balston Ultra High Purity (UHP) Nitrogen Generators are engineered to transform standard compressed air in to a safe regulated supply of 99.9999% pure nitrogen.

Typical applications include GC make up gas and carrier gas, including ECD (Electron Capture Detector), solvent evaporation, DSC (Differential Scanning Calorimeter) and virtually any analytical instrument that requires a small flow of ultra high puirty nitrogen.



Contact Information:

Parker Hannifin (UK) Limited **Process Air and Gas Business Unit**Hermitage Court, Hermitage Lane

Maidstone, Kent ME16 9NT

phone +44 (0)1622 723300 fax +44 (0)1622 728703 balstonukinfo@parker.com

www.parker.com/pag

Product Features:

- Produces a continuous supply of high purity nitrogen 99.9999% for GC applications
- Ideal for make-up and carrier gas applications including ECD
- Eliminate dangerous nitrogen cylinders from the laboratory
- Ideal for TOC applications including combustion, UV persulate and wet oxydation sampling techniques
- Low maintenance, minimal operation attention
- Designed to run 24 hours a day



ENGINEERING YOUR SUCCESS.

Nitrogen is produced by utilising a combination of filtration and pressure swing adsorption (PSA) technology. Standard compressed air is filtered by high efficiency coalescing filters to remove all contaminants down to 0.01 micron.

The air then passes through two columns filled with carbon molecular sieve (CMS) which adsorb O₂, CO₂, moisture and hydrocarbons. These are desorbed to atmosphere during the pressure swing cycle leaving a supply of ultra pure nitrogen.

For ultra sensitive applications such as ECD, units are also available with the addition of a heated catalyst module to ensure hydrocarbons are removed to < 0.1 ppm.

Principal Specification

Model	UHPN2-1100	HPN2-1100	HPN2-2000
Purity	99.9999 %	99.9999 %	99.99 %
Hydrocarbon Concentration	< 0.1ppm	N/A	N/A
Flow Rates	See Table	See Table	2,000 ml/min
Inlet Pressure	4.1 to 8.6 bar	4.1 to 8.6 bar	5.1 to 8.6 bar
Inlet Connection	1/4" NPT (Female)	1/4" NPT (Female)	1/4" NPT (Female)
Outlet Connection	1/8" NPT (Female)	1/8" NPT (Female)	1/8" NPT (Female)
Ambient Temperature	10 to 35°C	10 to 35°C	10 to 35°C
Electrical Requirements	230VAC - 50Hz	230VAC - 50Hz	230VAC - 50Hz
Power Consumption	700 Watts	25 Watts	25 Watts
Dimensions (H x W x D)	890 x 310 x 410 mm	890 x 310 x 410 mm	890 x 310 x 410 mm
Weight (Shipping)	47 Kg (50)	47 Kg (50)	47 Kg (50)

Flow Rates ml/min

Model	Inlet Air Pressure (Bar)	Maxiumum Outlet Flow (ml/min)	Maximum Outlet Pressure (Bar)
UHPN2-1100 and HPN2-1100	8.6	1,100	5.9
	7.6	1,000	5.2
	6.9	900	4.5
	6.2	800	4.1
	5.5	700	3.4
	4.8	600	3.1
	4.1	500	2.4

Ordering Information

Description	Model Number	
500 to 1,100 ml/min Ultra High Purity Nitrogen Generator (with Hydrocarbon removal)	UHPN2-1100	
500 to 1,100 ml/min High Purity Nitrogen Generator	HPN2-1100	
2,000 ml/min High Puirty Nitrogen Generator	HPN2-2000	
Installation Kit	IK7694	

Maintenance Items	Model Number	Change Frequency
Maintenance Kit for UHPN2-1100	MK7694	12 Months
Maintenance Kit for HPN2-1100	MK7692	12 Months
Maintenance Kit for HPN2-2000	MK7692	12 Months

© 2007 Parker Hannifin (UK) Limited

Distributed by:

\$3.2.026a January 2008





ADVANCED APPLIED TECHNOLOGIES

Contact Us:

Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com