

Dry Air Generators

for Nuclear Magnetic Resonance

Distributed by:



Carl Stuart Limited

ADVANCED APPLIED TECHNOLOGIES

Contact Us:

Irl Ph: 01 4523432

UK Ph: 08452 30 40 30

Web: www.carlstuart.com

Email: info@carlstuart.com

Ultra Dry Air
Gas Generator
Model UDA-300

NMR Gas on Demand, up to 340 lpm

The Parker Balston NMR (Nuclear Magnetic Resonance) gas generators are a complete system with carefully matched components engineered for easy installation, operation and long term reliability. They are designed to transform standard compressed air into a safe supply of dry (-73°C) gas ideal for ejecting, spinning and lifting on NMR instruments.

The NMR gas generators are an ideal alternative to nitrogen cylinders or dewars. Payback periods are typically less than one year.



Contact Information:

Parker Hannifin (UK) Limited
Process Air and Gas Business Unit
Hermitage Court, Hermitage Lane
Maidstone, Kent ME16 9NT

phone +44 (0)1622 723300
fax +44 (0)1622 728703
balstonukinfo@parker.com

www.parker.com/pag

Product Features:

- Produces a continuous supply of ultra dry (-73°C) air for NMR spectrometers and other analytical instruments
- Eliminate dangerous nitrogen cylinders or dewars from the laboratory
- Low cost alternative to nitrogen - pay back periods typically less than one year
- Compact, reliable and minimal maintenance
- Ideal gas supply for ejecting, spinning and lifting operations
- Designed to run 24 hours a day

 **BALSTON**
Analytical Gas Systems

ENGINEERING YOUR SUCCESS.

NMR gas is produced by utilising a combination of filtration and pressure swing adsorption (PSA) technologies. Standard compressed air is filtered by high efficiency coalescing filters to

remove all contaminants down to 0.01 micron. The air then passes through two columns filled with molecular sieve which adsorbs moisture. This is desorbed to atmosphere during the pressure

swing cycle leaving a supply of ultra dry NMR grade air. The PSA columns require no operator attention or maintenance. Simply connect the generator to the NMR instrument for consistent reliable analysis.

Principal Specification

Model	UDA-300
Dew Point	-73°C
Flow Rates	Up to 340 lpm
Inlet Pressure	4.1 to 8.6 bar
Inlet Connection	1/4" NPT (Female)
Outlet Connection	1/4" NPT (Female)
Ambient Temperature	10 to 35°C
Electrical Requirements	230VAC- 50Hz - 12VDC (to the UDA-300)
Power Consumption	10 Watts
Dimensions (H x W x D)	700 x 310 x 900 mm
Weight (Shipping)	93 Kg (98)

Ordering Information

Description	Model Number
NMR Gas Generator	UDA - 300
Installation Kit	IK7525

Maintenance Items	Model Number	Change Frequency
Maintenance Kit	MK7525	12 Months