


# Data sheet

## Resazurin cell viability assay

Cat. No: CA035  
10 000 assays (96-well format)

### Introduction

The **Resazurin cell viability assay** is a fluorescent assay that detects cellular metabolic activity. The kit offers a simple, rapid, reliable, sensitive, safe and cost-effective measurement of cell viability. Resazurin (**7-Hydroxy-3H-phenoxazin-3-one 10-oxide**) is a blue dye non-fluorescent until it is irreversibly reduced to the pink colored and highly red fluorescent resorufin by dehydrogenase enzymes in metabolically active cells. The fluorescent signal is monitored using 530-560 nm excitation wavelength and 590 nm emission wavelength. The absorbance is monitored at 570 nm and 600 nm. The fluorescent or colorimetric signal generated from the assay is proportional to the number of living cells in the sample.



### Kit Contents

Resazurin solution	4X25 ml
--------------------	---------

### Storage condition

Resazurin solution should be stored at  $-20^{\circ}\text{C}$  in the dark for long term storage. Thaw Resazurin solution completely and mix thoroughly before use. The product is stable for at least 10 freeze-thaw cycles. It can be stored at  $4^{\circ}\text{C}$  in the dark for up to 12 months.

### Assay procedure

1. Thaw out Resazurin solution (if kept frozen) and warm it to  $37^{\circ}\text{C}$  to ensure all components are completely in solution.
2. Plate cells into 96-well tissue culture plates using optimal cell concentration.
3. Carry out your experiment by adding agents of your interest into appropriate well and incubate with cells for a certain period of time.
4. Add Resazurin solution to plate (10% of the initial volume in the well). Return cells to the incubator and continue the incubation at for least 1 hour and up to 24 hours at  $37^{\circ}\text{C}$  (Incubation times may vary depending on the metabolic rates of the cell lines being tested).
5. Measure absorbance at 570 nm (If wavelength correction is available, set to 600 nm) **or fluorescence with excitation  $E_x=530-570$  nm and emission  $E_m=590-620$  nm** using a micro-titer plate reader.

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)

 **canvax**

## Notes

### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to [www.canvaxbiotech.com](http://www.canvaxbiotech.com) for Material Safety Data Sheet of the product.

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: [info@labunlimited.com](mailto:info@labunlimited.com)  
Web: [www.labunlimited.com](http://www.labunlimited.com)

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: [info@labunlimited.co.uk](mailto:info@labunlimited.co.uk)  
Web: [www.labunlimited.co.uk](http://www.labunlimited.co.uk)