Bedienungsanleitung / Operating Manual

Rührsystem / Stirring System

Distributed by:

IS 6 – Var 🤇

A DVANCED APPLIED TECHNOLOGIES Contact Us: Irl Ph: 01 4523432 UK Ph: 08452 30 40 30 Web: www.carlstuart.com Email: info@carlstuart.com

Carl Stuart Limited



BedienungsanleitungSeite 2Operating ManualPage 8

Contents

General information	8
Safety guidelines	9
Overview	10
Commissioning / Cleaning	11
Technical data	12
What to do if	13
Spare parts	13
Accessories	13

General information

Symbols used



indicates notes that you must read - for your own safety, the safety of others and to protect your meter from being damaged.



indicates notes that inform you of special features.

Accuracy when going to press

The use of advanced technology and the high quality standard of our instruments are the result of continuous development. This may result in differences between this operating manual and your instrument. We cannot guarantee that there are absolutely no errors in this manual. We are sure you will understand that we cannot accept any legal claims resulting from the data, figures or descriptions.

Warranty

The designated instrument is covered by a warranty of 12 months from the date of purchase.

The instrument warranty extends to manufacturing faults that are determined within the period of warranty. The warranty claim extends to restoring the instrument to readiness for use but not, however, to any further claim for damages. Improper handling or unauthorized opening of the instrument invalidates any warranty claim.

To ascertain the warranty liability, return the instrument and proof of purchase together with the date of purchase freight paid or prepaid.

Safety guidelines



Read these safety guidelines carefully before putting the instrument into operation!

This instrument has been built and tested according to the IEC 1010 safety standards for electronic measuring instruments and has left our works in a condition complying with all the requirements of technical safety.

The perfect functioning and operational safety of the instrument can only be ensured if the user observes the normal safety precautions as well as the specific safety guidelines stated in the present operating instruction.



ATTENTION: MAGNETISM! Effects of the magnetic field have to be taken into account (e.g. with cardiac pacemakers, data carriers...).

- Before connecting the plug-in power supply unit to the electricity supply network it must be ensured that the operating voltage stated on the plug-in power supply unit corresponds to the mains voltage (statement of the supply voltage range).
- The perfect functioning and operational safety of the instrument can only be maintained under the climatic conditions specified in the "Technical data" section of these operating instructions.
- When the instrument is moved from cold to warm surroundings, condensate may occur and interfere with the functioning of the instrument. In such a case, the user should wait until the temperature of the instrument has adapted to the ambient temperature before using the instrument again.
- Balancing, maintenance and repair work must only be carried out by a suitably qualified technician authorized by us.
- If there is reason to assume that the instrument can no longer be employed without a risk, it must be set aside and appropriately marked to prevent further use.
- The safety of the user may be endangered, e.g., if the instrument
 - shows visible damage,
 - no longer operates as specified,
 - has been stored over a longer period under unsuitable conditions,
 - has been subjected to difficult conditions during transport.
- If in doubt, the instrument should as a rule be sent back to the manufacturer's -"Wissenschaftlich-Technische Werkstätten GmbH" - for repair and maintenance.


Overview

The IS 6 – Var Stirring System has been especially developed for the WTW OxiTop[®]. All commercial sample bottles up to max. 93 mm bottom diameter can be used.

The stirring procedure is activated by a magnetic rotary field. This rotary field activates the permanently magnetic stirring rods.

The electronic interval control of the rotary field allows

- To place sample bottles on the stirring system at any time
- Disturbed stirring procedures to restart independently while the measuring operation is running.



- 1 Stirring power switch-over
- 2 Line adaptor connection

For BOD measurements with OxiTop[®] systems the stirring power should be set to "Low Power", as the higher power causes the temperature of the sample to rise.

Commissioning



Use the original line adaptor only.

The mains voltage must correspond to the operating voltage stated on the line adaptor!

- Place the stirring system into a thermostat cabinet or a thermostat box.
- Plug the line adaptor into the mains socket.
- Connect the line adaptor to the instrument (2).
 After a waiting time of approx. 3 s the inductive rotary field is automatically built up and thus the stirring rods are activated.



To guarantee an optimal and secure stirring function, use only original WTW stirring rods, RST 600.

- Place the prepared sample bottles on the center of the stirring places designed.
- Observe stirring operation, if necessary correct bottle position.

Operating note

1

Via the *Stirring power* switch-over (1) you can adjust the stirring power to meet the requirements of the sample.

When working with increased stirring power ("High Power"), the electronically controlled rotary field is operated with increased power. This causes an increased self-warming. If a perfect stirring operation is possible with reduced stirring power, work with "Low Power" preferably.

Cleaning

- Clean with a soft cloth and aqueous soap solution.



Do not use solvents or detergents that contain solvents (e. g. alcohol, acetone) for cleaning.

Technical data

Function principle	Magnetic rotary field		
Power supply (Line adaptor) <i>Euro version:</i> Friwo FW6798/11.8369, Friwo Part. No. 1769278 Input: 230 VAC (+10%/-15%)/ 50/60 Hz / 24 VA Output: 18 VAC / 1 A	USA version: Friwo FW6798/11.8371, Friwo Part. No. 1769294 Input: 120 VAC (+10%/-15%)/ 50/60 Hz / 23 VA Output: 18 VAC / 1 A		
Power consumption	24 VA max.		
Number of stirring places	6		
Diameter of vessel bottom	max. 93 mm		
Speed of stirring rods	Program-controlled 180 450 min ⁻¹		
Instrument security	Protection class: 3, IEC 1010		
	Protection degree: IP30		
EMC	Emissions: EN50081-1/FCC class A		
	Immunity: EN50082-1 EN50082-2, Namur		
Climate class	2, VDI/VDE 3540		
Ambient temperature	Storage: - 25 °C + 65 °C Operation: + 5 °C + 40 °C		
Relative humidity	Yearly mean: < 75 % 30 days/year: < 95 % Other days: < 85 %		
Test mark	CE		
Dimensions (WxDxH)	350 x 266 x 70 [mm]		
Weight (including line adaptor)	3.2 kg		

What to do if...

Error	Cause	Action		
No stirring function	 Mains power missing 	 Connect stirring system to power supply via line adaptor 		
	 Line adaptor defective 	 Replace line adaptor 		
	 Stirring system defective 	 Send stirring system to service department 		
Stirring function not sufficient	 Stirring rods unsuitable 	 Use original WTW stirring rods 		
	 Vessel bottom weakens the magnetic field (vessel bottom is too thick, crooked, or deformed) 	 Use suitable vessel 		
Centering of stirring rods faulty	 Vessel is not placed centrically on the stirring system 	 Adjust vessel position 		
	 Vessel bottom is too thick, crooked, or deformed 	 Use suitable vessel 		
	 Magnetic field is too weak 	 Send stirring system to service department 		

Spare parts

	Model	Order no.
Stirring rods	RST 600	209 120
Stirring rod remover	REF 600	209 130

Accessories

ĭ

The accessories that can be supplied are described in the WTW general catalog and in the price list.