

KU-3 Viscometer

for Paints, Coatings and Inks

ASTM D562 Compatible
(Industry specification)

Magnetic Spindle Coupling

Accommodates quart cans

Adapter provided
For pint and half pint cans

Easy to use
No weights, simplifies an
established test procedure

LED Display Info:

- Krebs Units
- Gram Units (Weight)
- Centipoise*

Lock-In Test Results
With hold key

Accuracy: $\pm 1.0\%$ of range

Repeatability: $\pm 0.5\%$

Standard Krebs Spindle

Measurement range:
40 to 141 KU, 32 to 1099 gm,
and 27 to 5274 cP*



What's Included

- Instrument
- Krebs-type Spindle (KU-1030)
- Adapter for Pint and Half-Pint Can

Note: Base holds Quart cans- No adapter needed

Optional Accessories

- Paste Spindle (KU-1031)

Applications

- Paints
- Coatings
- Adhesives
- Inks
- Pastes

*Centipoise values based on the conversion from Krebs Units as defined in ASTM D562.