

MULTIFARIOUS POSSIBILITIES: HEI-SHAKE SHAKERS

MADE IN
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

Multi Reax

Mixing powerfully – increasing sample throughput – saving time



The 3 mm shaking orbit in combination with a speed of 150 – 2,000 rpm guarantees the best mixing results

Two different adapters can be attached to the vortexer: for 2 tubes with 16 – 32 mm Ø or 6 tubes with 10 – 16 mm Ø

Comfortable inlets guarantee a secure hold when transporting the device

The speed can be adjusted via the speed controller. The parameters are shown precisely on the digital display

The timer allows the accomplishment of process steps efficiently

Distributed By:

Lab Unlimited
CARL STUART GROUP

Carl Stuart Limited
Tallaght Business Park
Whitestown
Dublin 24
D24 RFK3
Ireland

Ph: 00353 1 4523432
Fax: 00353 1 4523967
E-mail: onlinestore@labunlimited.com

Carl Stuart UK Limited
Calibre Scientific UK
Unit 5A, R-Evolution @
The Advanced Manufacturing Park
Selden Way
Rotherham
S60 5XA
United Kingdom

Ph: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: onlinestore@labunlimited.co.uk

heidolph
research made easy