




THE NEW EXHAUST FILTER

IMPROVEMENTS
ANYONE CAN SEE!


- Optimized design
- Consistently red
- Standard with splash protection and change indicator
- New development for vapor absorption and ventilation



**THE NEW EXHAUST FILTER.
AS OF NOW IN RED AND
STANDARD WITH SPLASH
PROTECTION AND
CHANGE INDICATOR.**



**ENLARGED INTAKE PORT.
IMPROVES THE ABSORPTION
OF SOLVENT VAPORS.**



**OPTIMIZED CASE DESIGN IMPROVES
AIR FLOW AND CIRCULATION OF
HAZARDOUS VAPORS THROUGH
OUR MULTICOMPONENT GRANULATE.**



**OPTIMIZED VENTILATION SLOTS
ENSURE CONSTANT FLOW
OF FILTERED EXHAUST AIR.**



EXHAUST FILTERS AS OF NOW CONSISTENTLY RED.

AS ALWAYS - PRICE ADVANTAGES WITH OUR USEFUL ECONOMY PACKS.



Part No.	Filter size	Recommended for waste containers with a volume of...	Features	Service life
610 534	S (small)	up to 5 liters	Splash protection with change indicator	3 Months
610 535	M (medium)	up to 20 liters	Splash protection with change indicator	6 Months
107 986	L (large)	more than 20 liters	Splash protection with change indicator	6 Months

Part No.	Description	Filter sizes	Number/unit
190 335	Economy package, exhaust filters with change indicator	S (small)	4 each
190 336	Economy package, exhaust filters with change indicator	M (medium)	3 each


Exhaust filters made from fire-resistant material and xl-versions for barrels you will find in our main catalogue and at www.scat-europe.com



ALL EXHAUST FILTERS YOU WILL FIND IN OUR NEW CATALOGUE, SAFETY SOLUTIONS 2016-2017



AND ONLINE AT WWW.SCAT-EUROPE.COM



facebook.com/SCATEurope



[@SCATEurope](https://twitter.com/SCATEurope)



SCAT Europe GmbH

www.scat-europe.com



Distributed by:

Lab Unlimited
CARL STUART GROUP

Tallaght Business Park
Whitestown, Dublin 24,
Ireland
D24 RFK3

Tel: (01) 4523432
Fax: (01) 4523967
E-mail: info@labunlimited.com
Web: www.labunlimited.com

Quatro House, Frimley Road,
Camberley,
United Kingdom
GU16 7ER

Tel: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: info@labunlimited.co.uk
Web: www.labunlimited.co.uk