

# Accessories for all Chambers and Glove Boxes: Stainless Steel Support Carts (Trolleys)



800-CART III

These carts are very useful for increasing portability and laboratory space.

All carts have anti-scoff casters including brakes and locks. Carts also feature a cushioning system that minimizes vibration, shocks, and wear and tear on chambers.

Most carts feature lower horizontal gas tank support racks which provide greater tank handling safety and flexibility.

Also available in a “knocked down” (unassembled) version. See information below.



800-CART/GB

PRODUCT NUMBER	CART DESCRIPTION	SIZES (O.D.)
CART-GB	For #GB, #NB, #AC Units	51”w x 24”d x 35” h
CART-GBB	For #GBB, #NBB, #ACB Units	55”w x 34”d x 35” h
CART-ABC	For #ABC Glove Box	43”w x 24”d x 34” h
CART-ABB	For #ABB Glove Box	49”w x 30”d x 34” h
CART-ABD	For #ABD Glove Box	63”w x 31”d x 34” h
CART-ABE	For #ABE Glove Box	76”w x 32”d x 34” h
CART-CLC	For #870 sized units	55”w x 34”d x 31” h

## Stainless Steel “Knocked Down” Cart (Trolley)

The Knocked Down Cart (Trolley) size is 40” (1016 mm) w x 24” (610 mm) d x 34” (864 mm) h.



Kit includes: Main base frame, casters, uprights, and top.



Insert casters and tighten



Force fit the four legs (uprights)



Place top on for legs & press into place.



800-CART KD

Distributed by:

**Lab Unlimited**  
CARL STUART GROUP

Tallaght Business Park  
Whitestown, Dublin 24,  
Ireland  
D24 RFK3

Tel: (01) 4523432  
Fax: (01) 4523967  
E-mail: info@labunlimited.com  
Web: www.labunlimited.com

Quatro House, Frimley Road,  
Camberley,  
United Kingdom  
GU16 7ER

Tel: 08452 30 40 30  
Fax: 08452 30 50 30  
E-mail: info@labunlimited.co.uk  
Web: www.labunlimited.co.uk