

Lab Unlimited
CARL STUART GROUP

A CALIBRE SCIENTIFIC COMPANY

Carl Stuart Limited

Tallaght Business Park
Whitestown
Dublin 24
D24 RFK3
Ireland

Ph: 00353 1 4523432
Fax: 00353 1 4523967
E-mail: onlinestore@labunlimited.com

Carl Stuart UK Limited

Calibre Scientific UK
Unit 5A, R-Evolution @
The Advanced Manufacturing Park
Selden Way
Rotherham
S60 5XA
United Kingdom

Ph: 08452 30 40 30
Fax: 08452 30 50 30
E-mail: onlinestore@labunlimited.co.uk

Julabo

THE TEMPERATURE CONTROL COMPANY

VALEGRO™

Recirculating Coolers

The new generation of chillers from JULABO
Compact – user-friendly – modern



**with natural
refrigerant**

ENGLISH



High level of user-friendliness and operational reliability

The ergonomically inclined OLED display and the easily accessible and illuminated level indicator ensure easy operation and enable placement under the lab bench. This saves time in the process and space in the work area. The removable venting grid warrants simple cleaning of the condenser and ensures long-term maintenance of performance and operational safety.

Intuitive menu navigation

All functions are simple and quick to access. Full-text instructions in German and English enable an immediate response to possible incidents and secure a safe handling of unit and application at all times. Useful functions such as timer, autostart and data recording improve process monitoring and lead to greater efficiency and time savings.



adjustable Pump

Adjustable pump pressure

VALEGRO recirculating coolers have an infinitely adjustable pump. The flow rate can therefore be optimally adapted to the application, which enables safe working even with pressure-sensitive glass apparatus and reduces the sound level at the same time.

Optimal connectivity

RS232 and USB-C interfaces are integrated as standard, located directly under the front display and easy to access. The optional Ethernet interface enables integration into any network and full access to all functions via a PC.



VALEGRO™ 350					VALEGRO™ 500			
BOersdteerl lN-Nor..					9610035			
AWrobrekitinstge tmeprmpreartautrubreer eriacn hg °eC °C					-20 ... +40			
Temperaturkeo sntsatsbainlizt y° C°C					±0..3			
KCäolotelilnegis ctuapngac °iCty °C					+20	0	-10	-20
(MBAetdhi uflmui dE tEhtahnaonlo) lk					0.18	0.25	0.18	0.07
Pumpe cnalpeiasctiutny gfl Foöwr dreartset rl/ommipl / min					19			
Pumpe cnalpeiasctiutny gp rFeörsduerred rbuacrkg					0.85			
Digitale i nStcehrnfaitctesstellen					RS232,, USB-C,, Ethernet (optionall)			
Menu functions					TimeTr-iFmuenrk-ftuionnc,t iPounm, lPnefinditreulcyk a sdtjufsetnalbolse epiumnstpe,lIbar, DRaetceonr da udfaztai,c hEtnheenr,n Eeth cearn ebte n activhtraätgeldig rhe tfroeiisspcehcatlitvb			
ATobtmale dssimunegnesnio gnes sWam×Dt ×BH×T c×mH cm					23 × 45 × 45			