

# LABORHITS

## Vortexer LLG-uniTEXER

**One head - multiples uses**

**Touch and permanent mixing modes**

**Three speeds levels**



### Features


- Touch function for short-term operation
- Carrier head is very easy to clean
- Precise speed control via a variable rotary switch
- Elastomeric feet ensure excellent stability on the bench and quiet operation
- Modern design
- Robust housing made of ABS



### Specifications

- Type of movement:** orbital/vortex  
**Orbital diameter:** 3,7 mm  
**Speed:** 1000, 2000 and 3000 rpm  
**Dimensions (W x D x H):** 173 x 198 x 198 mm  
**Weight:** 3,8 kg  
**Protection class:** IP 20  
**Power supply:** 220 V, 50/60 Hz  
**Warrenty:** 3 years

### 6 Different Sizes of Tubes, including standard assay plate



- a** Center pad area for hand held mixing
- b** 1 x Assay plate (to be used with retention cords)
- c** 2 x 50 ml conical tubes (horizontally mounted)
- d** 2 x 15 ml conical tubes (horizontally mounted)
- e** 2 x 5 ml tubes (horizontally mounted)
- f** in the same retention space as 15 ml) 4 x 1.5/2 ml microfuge tubes (vertically mounted)
- g** 6 x 0.5 ml microfuge tubes (vertically mounted)
- h** 24 x 0.2 ml microfuge tubes (strip tubes) (vertically mounted)

| Type         | PK | Cat. No.  |
|--------------|----|-----------|
| LLG-uniTEXER | 1  | 6.261 826 |